

Taking Cots to a New Level



The Impression Cot

A style of cot to meet the needs of customers in an ever more demanding environment. The cot is designed with ease of assembly and installation built in, enabling personal care and most importantly providing a safe environment of children of all ages and sizes.

Available in two mattress platform sizes:

1700x700mm External measurements 1940x840mm 1900x900mm External measurements 2140x1040mm

Cot surround and ends available at heights of: 1700x700mm - available in only 600 & 800mm 1900x900mm - available in 600,800 & 1000mm

Sides available in either Beech Bars or Plexiglas.

Full height padding in Oilcloth making cleaning easy and comes in a variety of colours.

The padding is attached to the cot frame with velcro making it easy to remove and clean.

Manufactured for ease of delivery and installation.

Technical Specifications		
Cot options Available:	1700x700x600H (mm)	1700x700x800H (mm)
	1900x900x600H (mm)	1900x900x800H (mm)
000° 20000° 20000° 20000° 20000° 20000° 20000° 2000	1900x900x1000H (mm)	
Weight:	Complete cot 80Kg	SWL 150Kg
Construction:	Cot framework platform	
NULL CONTRACTOR ACCOUNT ACCOUNT ACCOUNT ACCOUNT	& base	Steel st37
	Surface finish	Powder coated
Powered system	Linak CB-12 Motor height adjustable, Linak LA32. K-200	
Network	220-240V AC	

Feature

Electric height adjustment Bi folding doors on both sides

Electric backrest

Electric Trendelenburg & Anti Trendelenburg 40cm-80cm height adjustment

Safety latch on folding doors

Padding constructed in eight pieces

Strong steel construction 12.5 cm braking castor wheels

Knee Break (Optional)

Options Available



Beech Bars



Impression cot gate slider

Benefit

Improving moving and handling procedures.

Facilitates easy access to the child from either or both sides of the cot.

Facilitates easy raising of child into a sitting position, this assists in the reduction of respiratory problems. Allows tilt & anti tilt of cot if required.

Optimising the care environment whilst at the same time reducing the risk of back injuries. The adjustable height function is paramount in improving the delivery of therapy transfers and care.

The function can also be used to improve and maintain independence.

Helps maintain a safe environment preventing accidental or inappropriate opening.

Enabling easy delivery and easy removal of padding if necessary, also helps facilitate easy cleaning.

Ensures stability, durability and longevity.

Facilitates easy movement and cleaning whilst improving the safe environment of both child and caregiver. Facilitates full profiling of cot.



Impression cot showing standard and Vision Padding



SCAN MOBILITY LTD 222 Moss Ln, Hesketh Bank Lancashire PR4 6AE Tel: **01704 233344** Fax: 01704 233355 E-mail: sales@scanmobility.co.uk www.scanmobility.co.uk